

# 6inch Desk Model Aneroid Sphygmomanometer

**ACCOSON**

Established 1859  
Accuracy starts in the Name

The 6inch series is a range of aneroid sphygmomanometers offering the clinician the highest quality at an affordable price. A modern solution for all locations with innovative and practical features.

- 1** Made in UK and guaranteed accurate to ISO 81061-1
- 2** Range 0 – 300mmHg, with scale graduations from zero
- 3** New flexible shock mount resists damage to gauge
- 4** Clear white face with sharp black scale
- 5** Fine red pointer gives precise indication
- 6** Fitted with velcro style latex free cuff
- 7** Tough bright white ABS construction with built in carrying handle
- 8** Lightweight ABS cuff container, easy to clean

6inch Desk  
Code 0332



+44(0)1279 433456 | [accoson@accoson.com](mailto:accoson@accoson.com)

[www.accoson@accoson.com](http://www.accoson@accoson.com)

CE  
0473