

6inch Stand Model Aneroid Sphygmomanometer

ACCOSON

Established 1859
Accuracy starts in the Name

The 6inch series is a range of aneroid sphygmomanometers offering the clinician the highest quality at an affordable price. A modern solution for all locations with innovative and practical features.

- 1** Made in England and guaranteed accurate to ISO81060-1
- 2** Range 0 – 300 mmHg, with the scale calibrations starting at zero
- 3** Clear white face with fine red pointer gives pin point visibility
- 4** Smooth lines for easy cleaning without awkward areas
- 5** Designed to withstand double full scale overpressure
- 6** Adjustable height with inclined gauge
- 7** Stable base with five castors
- 8** Cuff stored behind the gauge in a large lightweight container with push handle
- 9** Latex free throughout with adult cuff and long coiled tubing



6inch Stand
Code 0362

+44(0)1279 433456 | accoson@accoson.com

www.accoson@accoson.com  0473

A C Cossor & Son (Surgical) Ltd.
Accoson Works, Parkway, Harlow Business Park, Harlow, Essex, CM19 5QP England

