


6inch Wall Model Aneroid Sphygmomanometer

ACCOSON

Established 1859
Accuracy starts in the Name

The 6 inch series is a range of aneroid sphygmomanometers offering the clinician the highest quality at an affordable price. A modern solution for all locations with innovative and practical features.

- 1 Made in England and guaranteed accurate to ISO 81060-1
- 2 Range 0 – 300 mmHg, with the scale calibrations starting at zero
- 3 Clear white face with fine pointer gives pin point visibility
- 4 Smooth lines for easy cleaning without awkward areas
- 5 Designed to withstand double full scale overpressure
- 6 Swivel wall mounting bracket
- 7 Long coiled tube supplied
- 8 Latex free throughout
- 9 Adult velcro cuff fitted



6inch Wall
Code 0342

+44(0)1279 433456 | accoson@accoson.com

www.accoson@accoson.com

CE
0473

A C Cossor & Son (Surgical) Ltd. Accoson Works, Parkway, Harlow Business Park, Harlow, Essex, CM19 5QP England