

Duplex Portable Hand Model Aneroid Sphygmomanometer

ACCOSON

Established 1859
Accuracy starts in the Name

1 Made in England and guaranteed accurate to ISO 81061-1

2 One handed use, with thumb and fingertip control

3 Lightweight diecast metal gauge, clear white face, robust and durable

4 Supplied in soft zipper case with coiled tube

5 Range 0 – 300mmHg with scale graduations from zero, withstands 100% overpressure

6 Fitted with ambidex anti-microbial cuff, with unique soft edges

Duplex Hand Model
Code 0322A



+44(0)1279 433456 | accoson@accoson.com

www.accoson@accoson.com

CE
0473

A C Cossor & Son (Surgical) Ltd. Accoson Works, Parkway, Harlow Business Park, Harlow, Essex, CM19 5QP England