

# Pocket Portable Aneroid Sphygmomanometer

ACCOSON

Established 1859  
Accuracy starts in the Name

**1** Made in UK and guaranteed accurate to ISO 81061-1

**2** Cuff/Pocket clip on the back of the lightweight diecast gauge

**3** Suitable for both left and right hands

**4** Supplied in soft zipper case

**5** Range 0 – 300mmHg with scale graduations from zero

**6** Clear white face with sharp black scale

**7** Fitted with ambidex anti-microbial cuff, with unique soft edges

Pocket Aneroid  
Code 0312A



+44(0)1279 433456 | [accoson@accoson.com](mailto:accoson@accoson.com)

[www.accoson@accoson.com](http://www.accoson@accoson.com)

CE  
0473