

Six⁰⁰ Series Desk Model Aneroid Sphygmomanometer

ACCOSON

Established 1859
Accuracy starts in the Name

1 Made in UK and guaranteed accurate to ISO 81061-1

2 Range 0 – 300mmHg, with the scale graduations starting at zero

3 Five year calibration warranty

4 New flexible shock mount resists damage to gauge

5 Clear white face with sharp black scale

6 Fine red pointer gives precise indication

7 Fitted with ambidex anti-microbial cuff, with unique soft edges

8 Tough bright white ABS construction with built in handle

9 Lightweight ABS cuff holder is easy to clean

Six00 Desk
Code 0632A



+44(0)1279 433456 | accoson@accoson.com

www.accoson@accoson.com **CE**
0473

A C Cossor & Son (Surgical) Ltd. Accoson Works, Parkway, Harlow Business Park, Harlow, Essex, CM19 5QP England