

# Six<sup>00</sup> Series

## Stand Model Aneroid Sphygmomanometer

**ACCOSON**

Established 1859  
Accuracy starts in the Name

The new SIX00 series is a premium range of aneroid sphygmomanometers offering the clinician the highest quality at an affordable price. A modern solution for all locations with innovative and practical features.

- 1** Made in England and guaranteed accurate to ISO 81060-1
- 2** Range 0 – 300 mmHg, with the scale calibrations starting at zero
- 3** Clear white face with fine red pointer gives pin point visibility
- 4** Attractive square profile with smooth lines for easy cleaning without awkward areas
- 5** Five year calibration warranty
- 6** Designed to withstand double full scale overpressure
- 7** Adjustable height with inclined gauge, and a stable base with five castors
- 8** Cuff stored behind the gauge in a large lightweight container with push handle
- 9** Latex free throughout with adult ambidex cuff and long coiled tubing



SIX00 Series  
Stand Model  
Code 0662A

+44(0)1279 433456 | [accoson@accoson.com](mailto:accoson@accoson.com)

[www.accoson@accoson.com](http://www.accoson@accoson.com)  0473

A C Cossor & Son (Surgical) Ltd.  
Accoson Works, Parkway, Harlow Business Park, Harlow, Essex, CM19 5QP England

